



---

- **PGP AI ML – Capstone Interim Submission**

- **Group A CV Batch 1**



Automating Pneumonia screening in chest radiographs and providing affected area details through powerful AI techniques can assist physicians to make better clinical decisions

- **Background** - Pneumonia accounts for over 15% of all deaths of children under 5 years old internationally. Accurately diagnosing pneumonia requires review of a chest radiograph (CXR) by highly trained specialists and confirmation through clinical history, vital signs and laboratory exams. Pneumonia usually manifests as an area or areas of increased opacity on CXR. However, the diagnosis of pneumonia on CXR is complicated because of a number of other conditions in the lungs such as fluid overload (pulmonary edema), bleeding, volume loss (atelectasis or collapse), lung cancer, or post-radiation or surgical changes. Outside of the lungs, fluid in the pleural space (pleural effusion) also appears as increased opacity on CXR. CXRs are the most commonly performed diagnostic imaging study. A number of factors such as positioning of the patient and depth of inspiration can alter the appearance of the CXR, complicating interpretation further. In addition, clinicians are faced with reading high volumes of images every shift. Tissues with sparse material, such as lungs which are full of air, do not absorb the X-rays and appear black in the image. Dense tissues such as bones absorb X-rays and appear white in the image
- **Problem statement** - We need to ***build a Pneumonia detection model*** to detect a visual signal for pneumonia in medical images. Specifically, the algorithm needs to automatically locate lung opacities on chest radiographs, providing affected area details through bounding box
- **Data Description:** Data provided is Dicom original images: - Medical images are stored in a special format called DICOM files (\*.dcm). They contain a combination of header metadata as well as underlying raw image arrays for pixel data. While we are theoretically detecting “lung opacities”, there are lung opacities that are not pneumonia related. In the data, some of these are labeled “Not Normal No Lung Opacity”. This extra third class indicates that while pneumonia was determined not to be present, there was nonetheless some type of abnormality on the image and oftentimes this finding may mimic the appearance of true pneumonia





# Exploratory Data Analysis

---








- Used pydicom library to read Dicom images
- 26684 unique patientIDs
- Few patients have multiple bounding boxes
  - 3266 patients have 2 bounding boxes defined
  - 119 patients have 3 bounding boxes defined
  - 13 patients have 4 bounding boxes defined
- We will run a binary classification to predict patients with pneumonia
- Target – 1 indicates patients with pneumonia and target – 0 indicates normal patients or patients having lung opacity but not pneumonia
- 8,851 (29.3%) patients are healthy/Normal
- 9,555 (31.6%) patients have Lung Opacity
- 11,821 (39.1%) patients have No Lung Opacity but may have other lung abnormality








	class	Target	Patient Count
0	Lung Opacity	1	9555
1	No Lung Opacity / Not Normal	0	11821
2	Normal	0	8851

- Patient distribution shows most pneumonia patients or patients are found between age 40-60
- More records for male is observed in the data
- Proportion of pneumonia and normal patients re equal for female patients
- However, for males, proportion of pneumonia patients are more than normal











**Lung opacity is seen  
maximum in age 50-65  
followed by age group 20-50**

- Chest X-ray can be done in **Posterior-Anterior (PA)** projection and **Anterior-Posterior (AP)** Projection
- Lung opacity seems to be more visible in AP position



## Visualization of sample X-ray images for each category

['Normal Xray | No Lung Opacity / Not Normal']




['Normal Xray | Normal']



['Normal Xray | No Lung Opacity / Not Normal']



['Pneumonia Infected | Lung Opacity']



['Normal Xray | No Lung Opacity / Not Normal']



['Normal Xray | No Lung Opacity / Not Normal']



['Pneumonia Infected | Lung Opacity']



['Pneumonia Infected | Lung Opacity']









['Normal Xray | No Lung Opacity / Not Normal']






Visualization of sample  
X-ray images with  
bounding boxes





# Model Building

---







Architectures	Accuracy	Precision	Recall	AUC-ROC
Basic CNN				
UNet				
Faster RCNN	78.4%	69%	15%	
VGG19	81%	55%	93%	91.7%
VGG16	80%	50%	81%	90.7%
ResNet50	77%	55%	80%	83.2%
InceptionNet v3	86%	55%	81%	86.8%

Since the problem statement requires us to read images, we use CNN and various state of the art deep learning architectures listed here to identify lung opacity with higher accuracy and recall

From the consolidated results, **VGG19** seems to be the best model with highest recall





Few cases of “Faster RCNN” model predicting bounding boxes of pneumonia patients with confidence intervals